

TM 2314 - SS SELECTIVE AGAR, IMPROVED, ((DOUBLE PACK))

INTENDED USE

Used for the selective isolation and differentiation of *Salmonella* and *Shigella* species.

PRODUCT SUMMARY AND EXPLANATION

SS Selective Agar, Improved is recommended as selective medium for the isolation of *Salmonella* as well a *Shigella* species from clinical specimens. It provides significantly greater sensitivity and specificity in the detection of both the organisms. The other selective medias like HE, SS and XLD largely fail to suppress the growth of *Salmonella* interfering organism like *Citrobacter* and *Proteus* which resemble the presence of *Salmonella*.

Bacteria that decarboxylate lysine to cadaverine can be recognized by the appearance of a red colouration around the colonies due to an increase in pH. Sodium thiosulphate is reduced by certain species of enteric organisms to sulphite and H₂S gas and this reductive enzyme process is attributed by thiosulphate reductase. Production of H₂S gas is detected as an insoluble black precipitate of ferrous sulphide, formed upon reaction of H₂S with ferric ions or ferric citrate, indicated in the centre of the colonies. Part B addition to the medium helps in improving the selectivity of the medium. It selectively inhibits the growth of gram positive organisms. The growth of *Proteus* species is also reduced.

COMPOSITION

Ingredients	Gms / Ltr
Part I	
Proprietary	81.93
Part II	
Proprietary	4.600

PRINCIPLE

The sugars xylose, lactose and sucrose provide sources of fermentable carbohydrates, xylose is mainly incorporated into the medium since it is not fermented by *Shigella* but practically by all enterics. This helps in the differentiation of *Shigella* species. Degradation of xylose, lactose and sucrose to acid causes phenol red indicator to change its colour to yellow. Addition of Cellobiose controls the growth of false positive *Salmonella* suspect. Lysine is included to differentiate the *Salmonella* group from the non-pathogens. *Salmonella* rapidly ferment xylose and exhaust the supply. Subsequently lysine is decarboxylated by the enzyme lysine decarboxylase to form amines with reversion to an alkaline pH that mimics the *Shigella* reaction.

INSTRUCTION FOR USE

- Dissolve 81.93 grams in 1000 ml distilled water. Add 4.6 ml of Part II.
- Heat to boiling with frequent agitation to dissolve the medium completely, do not autoclave or overheat. Overheating may destroy the selectivity of the medium.
- Cool to 45-50°C. Mix and pour into sterile Petri plates.

QUALITY CONTROL SPECIFICATIONS

Appearance of Powder	: Light yellow to pink homogeneous free flowing powder.
Appearance of prepared medium	: Reddish orange coloured clear to slightly opalescent gel forms in Petri plates.
pH (at 25°C)	: 7.4±0.2

INTERPRETATION



Cultural characteristics observed after incubation. Recovery rate is considered as 100% for bacteria growth on Soyabean Casein Digest Agar.

Microorganism	ATCC	Inoculum (CFU/ml)	Growth	Recovery	Color of the colony	Incubation Temperature	Incubation Period
<i>Salmonella</i> Typhimurium	14028	50-100	Luxuriant	≥70%	Red with black centers	35-37°C	18-48 Hours
<i>Escherichia coli</i>	25922	50-100	Fair-good	20-40%	Yellow	35-37°C	18-48 Hours
<i>Salmonella</i> Enteritidis	13076	50-100	Good-luxuriant	≥50%	Red with black centers	35-37°C	18-48 Hours
<i>Salmonella</i> Typhi	6539	50-100	Good-luxuriant	≥50%	Red with black centers	35-37°C	18-48 Hours
<i>Shigella dysenteriae</i>	13313	50-100	Good-luxuriant	≥50%	Red	35-37°C	18-48 Hours
<i>Shigella flexneri</i>	12002	50-100	Fair-good	20-40 %	Red	35-37°C	18-48 Hours
<i>Shigella sonnei</i>	25931	50-100	Fair-good	20-40 %	Red	35-37°C	18-48 Hours
<i>Klebsiella aerogenes</i>	13048	50-100	Fair	20-30 %	Yellow	35-37°C	18-48 Hours
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	25923	≥10 ³	Inhibited	0%	-	35-37°C	18-48 Hours
<i>Enterococcus faecalis</i>	29212	≥10 ³	Inhibited	0%	-	35-37°C	18-48 Hours

<i>Proteus mirabilis</i>	25933	50-100	None-poor	0-10%	Red	35-37°C	18-48 Hours
--------------------------	-------	--------	-----------	-------	-----	---------	-------------

PACKAGING:

In pack size of 500 gm bottles.

STORAGE

Dehydrated powder, hygroscopic in nature, store in a dry place, in tightly-sealed containers between 25-30°C and protect from direct sunlight. Under optimal conditions, the medium has a shelf life of 4 years. When the container is opened for the first time, note the time and date on the label space provided on the container. After the desired amount of medium has been taken out replace the cap tightly to protect from hydration.


Product Deterioration: Do not use if they show evidence of microbial contamination, discoloration, drying or any other signs of deterioration.

DISPOSAL

After use, prepared plates, specimen/sample containers and other contaminated materials must be sterilized before discarding.

REFERENCES

1. Isenberg, H.D. Clinical Microbiology Procedures Handbook 2nd Edition.
2. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.
3. Mallinson, E.T.et al.2000.Improved plating media for the detection of *Salmonella* species with typical and atypical hydrogen sulfide production.J.Vet. Diagn.Invest.12:83-87.
4. Mallinson, E.T.1991. Novelsal monella detection systemdeveloped; combines increased reliability, practicality. Feedstuffs 63:40-44.
5. Miller,R.G.et al.1991.Xylose-Lysine-tergitol 4:an improved selective agar medium for the isolation of *Salmonella*
6. Poult.Sci.70:2429-2432.Erratum, Poult.Sci.71:398,1992.
7. Miki,K.et al.1996.Re-speciation of the original strains of serovars in the *Citrobacter freundii* (Bethesda-Ballerup group) antigenic scheme of West and Edwards.Microbiol. Immunol.40:915-921.
8. Pollock, H.M. and B.J. Dahlgren.1974, Clinical evaluation of enteric media in the primary isolation of *Salmonella* and
9. *Shigella*. Appl.Microbiol.27:197-201.
10. Salfinger Y., and Tortorello M.L. Fifth (Ed.), 2015, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.



NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

***For Lab Use Only**
Revision: 08 Nov., 2019