

## TS 012 – BORDETELLA SELECTIVE SUPPLEMENT

### INTENDED USE

For selective isolation of *Bordetella pertussis*.

### COMPOSITION

Ingredients	Concentration
Cephalexin	20mg

(Per vial sufficient for 1000 ml medium)

### Applied in Medium

TM 044 - BORDET GENGOU AGAR BASE

TM 1147 - BORDET GENGOU AGAR BASE W/ 1.6% AGAR

TM 062 - CHARCOAL AGAR BASE













TM 063 - CHARCOAL AGAR BASE W/ NIACIN

### INSTRUCTION FOR USE

Rehydrate the content of 1 vial aseptically with 2 ml of sterile distilled water. Mix well and aseptically add it to 500 ml of sterile, molten, cooled (45-50°C) B TM 044 - BORDET GENGOU AGAR BASE / TM 1147 - BORDET GENGOU AGAR BASE W/ 1.6% AGAR or 1000 ml of sterile, molten TM 063 - CHARCOAL AGAR BASE W/ NIACIN together with 10% v/v defibrinated horse blood. Mix well and pour into sterile petri plates. The vial content may be added to 500 ml of sterile half strength TM 062 - CHARCOAL AGAR BASE with 10% v/v defibrinated horse blood for use as a transport medium for *Bordetella pertussis*.

### STORAGE

Vials should be stored in sealed container at 2-8°C.

 Good Manufacturing Practices Certified	 For In Vitro Diagnostic Use	 Quantity	 Lot / Batch Number	 Catalogue Number	 Manufacturer
 Temperature Unit	 Authorized Representative MedWor GmbH Buckhorn 10, 40143 Essen, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

**NOTE:** Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

**\*For Lab Use Only**  
Revision: 5 March., 2022