

TS 073 – PARK AND SANDERS SELECTIVE SUPPLEMENT A

INTENDED USE

For selective isolation and enumeration of thermotolerant *Campylobacter* species.

COMPOSITION

Ingredients	Concentration
Vancomycin	10mg
Trimethoprim	10mg

(per vial sufficient for 950 ml medium)

Applied in Medium

TM 1335 - CAMPYLOBACTER CEFEX BROTH BASE

TM 1403 - PARK & SANDER ENRICHMENT BROTH BASE













TM 2274 - PARK & SANDER ENRICHMENT BROTH BASE (ISO/DIS 10272:1995)

INSTRUCTION FOR USE

Rehydrate the contents of 1 vial aseptically with 5 ml of 50% ethanol. Mix well and aseptically add to 940 ml of sterile, cooled TM 1335 - CAMPYLOBACTER CEFEX BROTH BASE / TM 1403 - PARK & SANDER ENRICHMENT BROTH BASE / TM 2274 - PARK & SANDER ENRICHMENT BROTH BASE (ISO/DIS 10272:1995) along with 50 ml sterile defibrinated blood and Park and Sanders Selective Supplement B. Mix well and dispense into sterile tubes.

STORAGE

Vials should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 IVD For In Vitro Diagnostic Use	 QTY. Quantity	 LOT/ B. NO. Lot / Batch Number	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 EC REP Authorized Representative MedStar GmbH Buckhorn 10, 40143 Alsdorf, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

***For Lab Use Only**
Revision: 4 March., 2022