

TS 176 – RPF SUPPLEMENT (Fibrinogen plasma trypsin inhibitor Supplement)

INTENDED USE

For isolation and enumeration of coagulase positive staphylococci from food and clinical samples.

COMPOSITION

Ingredients	Concentration
Bovine fibrinogen	0.375g
Rabbit plasma	2.500ml
Trypsin inhibitor	2.500mg
Potassium tellurite	2.500mg

(per vial sufficient for 100 ml medium)

Applied in Medium

TM 358 - BAIRD PARKER AGAR BASE

TMV 358 - BAIRD PARKER AGAR BASE(VEG.)

TM 1579 - BAIRD PARKER AGAR BASE (RPF) (ISO 6888-1 & 2:1999)

TM 2300 - RPF AGAR BASE (ISO 6888-2:1999)

INSTRUCTION FOR USE













Rehydrate the contents of 1 vial by adding 10 ml of sterile distilled water. Turn the vial upside-down slowly two to three times avoiding frothing. Alternately the contents may be dissolved by gentle stirring for 1-2 minutes with a sterile glass rod or a stirrer. Immediate dissolution is not obtained. Keep it standing for 1-2 hours for better dissolution. Slight particulate matter may be observed which does not affect the performance parameter of the supplement. Aseptically add the supplement in 90 ml of sterile, molten TM 358 - BAIRD PARKER AGAR BASE / TMV 358 - BAIRD PARKER AGAR BASE(VEG.) / TM 1579 - BAIRD PARKER AGAR BASE (RPF) (ISO 6888-1 & 2:1999) / TM 2300 - RPF AGAR BASE (ISO 6888-2:1999) cooled to 48°C. Add the supplement slowly from the side of the flask with simultaneous gentle swirling of the flask to mix the supplement and the medium homogeneously. Pour into sterile petri plates and use immediately.

Note : The temperature of the molten, sterile Baird Parker Agar Base or Baird Parker Agar Base (Veg) should be strictly 48°C at the time of addition of the supplement. Any variation in this temperature may give erroneous false negative results.

STORAGE

Vials should be stored in sealed container at 2-8°C.



 GMP Good Manufacturing Practices Certified	 IVD For In Vitro Diagnostic Use	 QTY. Quantity	 LOT/ B. NO. Lot / Batch Number	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 EC REP Authorized Representative <small>MedWet GmbH Südstraße 10 48163 Aachen, Germany</small>	 European Conformity	 QR Code	 Consults instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

***For Lab Use Only**
Revision: 3 March., 2022