

TS 249 – MYCOPLASMA UROGENITAL SELECTIVE SUPPLEMENT

INTENDED USE

For selective cultivation of *Mycoplasma* associated with urogenital infections.

COMPOSITION

Ingredients	Concentration
Penicillin	100000Unit
Nalidixic acid	5mg
Amphotericin B	1mg

(Per vial sufficient for 500 ml medium)

Applied in Medium













TM 1267 - MYCOPLASMA UROGENITAL MEDIUM BASE

INSTRUCTION FOR USE

Aseptically rehydrate the contents of 1 vial in 10 ml sterile distilled water. Aseptically add 1 vial to 500 ml of sterile, cooled (45-50°C) TM 1267 - MYCOPLASMA UROGENITAL MEDIUM BASE along with 1 vial of Vitamins amino Growth Supplement, 1 vial of Urea Solution 5% and 50 ml Horse Serum. Mix well and dispense as desired.

STORAGE

Vials should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 IVD For In Vitro Diagnostic Use	 QTY. Quantity	 LOT/ B. NO. Lot / Batch Number	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 EC REP Authorized Representative MedWor GmbH Buckhorn 10, 40112 Düsseldorf, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

***For Lab Use Only**
Revision: 7 March., 2022