

TS 278 – CAMPYLOBACTER SELECTIVE SUPPLEMENT, ABEYTA

INTENDED USE

Recommended by FDA BAM for the selective isolation of *Campylobacter* species.

COMPOSITION

Ingredients	Concentration
Sodium cefoperazone	16.00mg
Rifampicin	5.00mg
Amphotericin B	1.00mg

(per vial sufficient for 500 ml medium)

Applied in Medium









TM 2105 - HEART INFUSION YEAST EXTRACT AGAR BASE

INSTRUCTION FOR USE

Rehydrate the contents of 1 vial aseptically with 1 ml of 0.2 N Sodium hydroxide, further add 4 ml of sterile distilled water. Mix well and aseptically add to 470 ml of sterile, molten, cooled (45-50°C) TM 2105 - Heart Infusion Yeast Extract Agar Base along with 25 ml of sterile lysed horse blood and reconstituted contents of one vial of Campylobacter Growth Supplement. Mix well and pour into sterile Petri plates.

STORAGE

Vials should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 IVD For In Vitro Diagnostic Use	 QTY. Quantity	 LOT/ B. NO. Lot / Batch Number	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 EC REP Authorized Representative MedNet GmbH Buckhorn 16, 40143 Düsseldorf, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

***For Lab Use Only**

Revision: 11 Feb., 2022