

## TBL 056 - STERILE DISCS

### INTENDED USE

Sterile discs can be used for testing antibiotic susceptibility testing by loading antibiotic of choice.

### QUALITY CONTROL SPECIFICATIONS

**Appearance** : Sterile filter paper discs of 6mm diameter.

### INTERPRETATION

Sterility of the discs checked for 14 days in Alternative Thioglycollate Medium (M010) at 30-35°C and Soyabean Casein Digest Medium (TMTH 011) at 20-25°C.

Sample	Results in Alternative Thioglycollate Medium	Results in Soyabean Casein Digest Medium
Test	Clear, No growth observed	Clear, No growth observed
Negative Control	Clear, No growth observed	Clear, No growth observed
Positive control	Luxuriant growth observed	Luxuriant growth observed

### PACKAGING:

In pack size of 50 Discs/vl.

### STORAGE

Store at 2 - 8°C. Use before expiry date on the label.

 Good Manufacturing Practices Certified	 For In Vitro Diagnostic Use	 Quantity	 Lot / Batch Number	 Catalogue Number	 Manufacturer
 Temperature Unit	 Authorized Representative	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

**NOTE:** Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

\*For Lab Use Only  
Revision: 08 Nov., 2019

