

TK 005– ALBERT'S METACHROMATIC STAIN – KIT

REF NO: TBL/QA/SS/LC/TK 005	STANDARD SPECIFICATION FOR ALBERT'S METACHROMATIC STAIN - KIT	Rev. No: 04
Issue No: 05		Review Date: 02/02/2028
Issue Date: 03/02/2025		Product code: TK 005

PRODUCT PROPERTIES

C.A.S Number	NA
Chemical Formula	NA
Formula weight	NA
Functional Uses	Albert's Metachromatic Stain - Kit is used for observing metachromatic granules
Standard Packaging	1 Kit (125 ml)
Key Component	<ul style="list-style-type: none"> Albert's Stain A Albert's Stain B

PHYSICAL PARAMETERS

PRODUCT PARAMETER	SPECIFICATION
Sensitivity test	This product has been tested and conforms to quality standard.

INSTRUCTION FOR USE

STAINING PROCEDURE
<ul style="list-style-type: none"> Prepare a thin smear on a clean dry glass slide.. Allow it to dry and fix with gentle heat. Stain with Alberts Stain A for 3-5 minutes. Drain the solution, do not wash. Apply Alberts Stain B for 1 minute. Rinse with water, blot dry. Examine the smear under microscope, 40X and 100X in oil immersion.

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

