

TM 1213 – LPM AGAR BASE

INTENDED USE

For isolation and cultivation of *Listeria monocytogenes* from food and dairy products.

PRODUCT SUMMARY AND EXPLANATION

L. monocytogenes is a gram-positive foodborne human pathogen responsible for serious infections in pregnant women that may ultimately result in abortion, stillbirth, birth of child with meningitis or primary bacteremia in adults and juveniles. The organism has been isolated from commercial dairy and other food processing plants, and is ubiquitous in nature, being present in a wide range of unprocessed foods and in soil, sewage and silage and river water. *Listeria* species grow over a pH range of 4.4 to 9.6, and survive in food products with pH levels outside these parameters. Motility is most pronounced at 20°C.

Lee and McClain developed LPM Agar, which is a modification of McBride *Listeria* Agar. It enhances the recovery of low numbers of *L. monocytogenes* from mixed microflora in samples. APHA also recommends this medium for food and dairy sample testing.

COMPOSITION

Ingredients	Gms / Ltr
Tryptone	5.000
Peptone	5.000
Beef extract	3.000
Glycine anhydride	10.000
Lithium chloride	5.000
Sodium chloride	5.000
Phenylethyl alcohol	2.500
Agar	15.000

PRINCIPLE

This medium consists of Tryptone, Peptone and Beef extract which are sources of nitrogen, vitamins and minerals. Sodium chloride maintains the osmotic balance of the medium. Glycine anhydride improves recovery of *Listeria*. Lithium chloride, moxalactam and phenyl ethanol aids in suppression of both gram-positive and gram-negative organisms including *Staphylococcus*, *Proteus* and *Pseudomonas* species. *Listeria monocytogenes* show blue-green iridescence when examined with oblique transmitted light.

INSTRUCTION FOR USE

- Dissolve 50.50 grams in 1000 ml purified/distilled water.
- Heat to boiling to dissolve the medium completely.
- Sterilize by autoclaving at 15 psi pressure (121°C) for 15 minutes.
- Cool to 45-50°C and aseptically add rehydrated contents of 1 vial of Moxalactam Supplement.
- Mix well and pour into sterile Petri plates.

QUALITY CONTROL SPECIFICATIONS



Appearance of Powder : Cream to light yellow homogeneous free flowing powder.
Appearance of prepared medium : Light yellow coloured clear to slightly opalescent gel forms in Petri plates.
pH (at 25°C) : 7.3 ± 0.2

INTERPRETATION

Cultural characteristics observed with added Moxalactum Supplement, after incubation.

Microorganism	ATCC	Inoculum (CFU/ml)	Growth	Recovery	Incubation Temperature	Incubation Period
<i>Escherichia coli</i>	25922	$\geq 10^3$	Inhibited	0%	35-37 °C	24 - 48 Hours
<i>Listeria monocytogenes</i> serovar 1	19111	50-100	Good-luxuriant	$\geq 50\%$	35-37 °C	24 - 48 Hours
<i>Listeria monocytogenes</i>	19112	50-100	Good-luxuriant	$\geq 50\%$	35-37 °C	24 - 48 Hours
<i>Listeria monocytogenes</i>	19117	50-100	Good-luxuriant	$\geq 50\%$	35-37 °C	24 - 48 Hours
<i>Pseudomonas aeruginosa</i>	27853	$\geq 10^3$	Inhibited	0%	35-37 °C	24 - 48 Hours
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	25923	$\geq 10^3$	Inhibited	0%	35-37 °C	24 - 48 Hours

PACKAGING:

In pack size of 100 gm bottles.

STORAGE

Dehydrated powder, hygroscopic in nature, store in a dry place, in tightly-sealed containers between 25-30°C and protect from direct sunlight. Under optimal conditions, the medium has a shelf life of 4 years. When the container is opened for the first time, note the time and date on the label space provided on the container. After the desired amount of medium has been taken out replace the cap tightly to protect from hydration.

Product Deterioration: Do not use if they show evidence of microbial contamination, discoloration, drying or any other signs of deterioration.

DISPOSAL













After use, prepared plates, specimen/sample containers and other contaminated materials must be sterilized before discarding.

REFERENCES

1. Beams R. E. and Girard K. F., 1959, Am. J. Med. Technol., 25:120.
2. Bremer P. and Osborne C., 1995, J. Food Prot., 58:604.
3. Isenberg, H.D. Clinical Microbiology Procedures Handbook 2nd Edition.
4. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.



5. Lee W. H. and McClain D., 1986, Appl. Environ. Microbiol., 52:1215.
6. Marshall R. T., (Ed.), 1992, Standard Methods for the Examination of Dairy Products, 16th Ed., APHA, Washington, D.C.
7. Murray P. R., Baron J. H., Pfaller M. A., Tenover J. C. and Tenover J. C., (Ed.), 2003, Manual of Clinical Microbiology, 8th Ed., American Society for Microbiology, Washington, D.C.
8. Patel J. R., Hwang C. A., Beuchat L. R., Doyle M. P. and Brackett R. E., 1995, J. Food Prot., 58:244
9. Salfinger Y., and Tortorello M.L., 2015, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.
10. Wehr H. M., 1987, J. Assoc. Off. Anal. Chem., 70:769.

 GMP Good Manufacturing Practices Certified	 IVD For In Vitro Diagnostic Use	 QTY. Quantity	 LOT/ B. NO. Lot / Batch Number	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 EC REP Authorized Representative <small>MedNet GmbH Birkstrasse 10 48163 Münster, Germany</small>	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

***For Lab Use Only**
Revision: 08 Nov., 2019