

TMC 005 - STERILE DISPOSABLE LOOP (P.P)

INTENDED USE

Recommended for to pick up and transfer a small sample (inoculum) from a culture of microorganisms..

Parameters	Size
Inner Diameter	3.5 mm
Out Diameter	6.0 mm
Loop calibrated	10µl

APPLICATIONS:


- Tools to safely transfer microbial specimen.
- Used in Clinical, Pharmaceuticals, testing Labs.

QUALITIES:

- The Loop do not requires flaming for sterilization.
- Loops are made from a plastic resin that provides controlled flexibility..
- They resist cutting the agar surface yet provide sufficient firmness to pick stubborn colonies. Loops are color-coded for quick ID.

PRODUCT FEATURES:

- Sterile
- Pack Size: 5 x 10 Nos & 1 x50 Nos.



NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

***For Lab Use Only**

