

TMC 050 - ASSORTED NICHROME LOOP

INTENDED USE

An inoculation loop, simple tool used for streaking. The loop of wire at the tip is made of nichrome and handle with stainless steel.

Parameters	Size
Diameter	4 mm, 2.5 mm, 1.3 mm
Calibrated to	0.01 ml, 0.05 ml, 1 µl

PRODUCT FEATURE:


Heat Resistant handle.

APPLICATION:

Microbiology laboratory, Molecular biology laboratory, Culture preparation and various other laboratories.

PACKAGING:

- Pack Size: 50 Pcs.



NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

***For Lab Use Only**
Revision: 23rd Jan 2023

