

TS 016 – LEGIONELLA SELECTIVE SUPPLEMENT-II

INTENDED USE

For selective isolation of *Legionella* species from mixed cultures.

COMPOSITION

Ingredients	Concentration
Cephalothin	2mg
Colistin sulphate	8mg
Vancomycin	0.250mg
Cycloheximide	40mg

(Per vial sufficient for 500 ml medium)

Applied in Medium

TM 684 - BUFFERED CHARCOAL YEAST EXTRACT AGAR BASE

INSTRUCTION FOR USE

Rehydrate the contents of 1 vial aseptically with 5 ml of 50% ethanol. Mix well and aseptically add to 500 ml of sterile, molten, cooled (45-50°C) TM 684 - BUFFERED CHARCOAL YEAST EXTRACT AGAR BASE along with rehydrated contents of Legionella Selective Supplement IV and Legionella Supplement TS 019. Mix well and pour into sterile petri plates.

STORAGE

Vials should be stored in sealed container at 2-8°C.

 Good Manufacturing Practices Certified	 For In Vitro Diagnostic Use	 Quantity	 Lot / Batch Number	 Catalogue Number	 Manufacturer
 Temperature Unit	 Authorized Representative MedNet GmbH Buckhorn 10, 48147 Münster, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

***For Lab Use Only**

Revision: 7 March., 2022