

TS 055 – NON SPORE ANAEROBIC SUPPLEMENT

INTENDED USE

For selective isolation or cultivation of non-sporing anaerobic bacteria.

COMPOSITION

Ingredients	Concentration
Sodium pyruvate	500mg
Menadione	0.250mg
Ferric pyrophosphate, soluble	2.500mg
Nalidixic acid	5mg

(per vial sufficient for 500 ml medium)

Applied in Medium

TM 915 - WILKINS CHALGREN ANAEROBIC AGAR BASE

TM 916 - WILKINS CHALGREN ANAEROBIC BROTH BASE

INSTRUCTION FOR USE

Rehydrate the contents of 1 vial aseptically with 10 ml sterile distilled water. Mix well and aseptically add to sterile, molten, cooled (45-50°C) 490 ml of TM 915 - WILKINS CHALGREN ANAEROBIC AGAR BASE / TM 916 - WILKINS CHALGREN ANAEROBIC BROTH BASE, enriched with 5% v/v sterile defibrinated blood. Mix well and pour into sterile petri plates.

STORAGE

Vials should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 IVD For In Vitro Diagnostic Use	 QTY. Quantity	 LOT/ B. NO. Lot / Batch Number	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 EC REP Authorized Representative MediStar GmbH Buckhorn 16, 48143 Münster, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

***For Lab Use Only**
Revision: 4 March., 2022