

TS 056 – G.N. SPORE ANAEROBIC SUPPLEMENT

INTENDED USE

For the selective isolation of gram-negative anaerobes.

COMPOSITION

Ingredients	Concentration
Sodium succinate	1250mg
Hemin	2.500mg
Vancomycin	1.250mg
Menadione	0.250mg
Nalidixic acid	5mg

(per vial sufficient for 500 ml medium)

Applied in Medium

TM 915 - WILKINS CHALGREN ANAEROBIC AGAR BASE

TM 916 - WILKINS CHALGREN ANAEROBIC BROTH BASE

INSTRUCTION FOR USE

Rehydrate the contents of 1 vial aseptically with 10 ml sterile distilled water. Mix well and aseptically add it to sterile, molten, cooled (45-50°C) 490 ml of TM 915 - WILKINS CHALGREN ANAEROBIC AGAR BASE / TM 916 - WILKINS CHALGREN ANAEROBIC BROTH BASE enriched with 5% v/v sterile defibrinated blood. Mix well and pour into sterile petri plates.

STORAGE

Vials should be stored in sealed container at 2-8°C.

 Good Manufacturing Practices Certified	 For In Vitro Diagnostic Use	 Quantity	 Lot / Batch Number	 Catalogue Number	 Manufacturer
 Temperature Unit	 Authorized Representative MedWise GmbH Buckhorn 10 40142 Düsseldorf, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

***For Lab Use Only**
Revision: 28 Feb., 2022