

TS 115 – LEGIONELLA (GVPC) SELECTIVE SUPPLEMENT

INTENDED USE

For isolation of *Legionella* species.

COMPOSITION

Ingredients	Concentration
Glycine	1.500g
Vancomycin hydrochloride	0.500mg
Polymyxin B sulphate	40000IU
Cycloheximide	40mg

(Per vial sufficient for 500 ml medium)

Applied in Medium

TM 1218 - LEGIONELLA AGAR BASE

TM 1894 - BUFFERED CHARCOAL YEAST EXTRACT AGAR MEDIUM (ISO 11731-2:2017)

TM 2156 - LEGIONELLA AGAR BASE W/O CHARCOAL

INSTRUCTION FOR USE

Rehydrate the contents of 1 vial aseptically with 10 ml of sterile distilled water. Mix well and aseptically add to 440 ml sterile, molten, cooled (45-50°C) TM 1218 - LEGIONELLA AGAR BASE / TM 2156 - LEGIONELLA AGAR BASE W/O CHARCOAL along with 50 ml of Legionella Growth Supplement BCYE or to 500 ml sterile, molten, cooled (45-50°C) TM 1894 - BUFFERED CHARCOAL YEAST EXTRACT AGAR MEDIUM (ISO 11731-2:2017). Mix well and pour into sterile petri plates.

STORAGE

Vials should be stored in sealed container at 2-8°C.

 Good Manufacturing Practices Certified	 For In Vitro Diagnostic Use	 Quantity	 Lot / Batch Number	 Catalogue Number	 Manufacturer
 Temperature Unit	 Authorized Representative MedWet GmbH Buckhorn 16, 48143 Münster, Germany	 European Conformity	 QR Code	 Consults instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

***For Lab Use Only**
Revision: 7 March., 2022