

TS 171 – RIPPEY CABELLI SUPPLEMENT

INTENDED USE

For cultivation of *Aeromonas hydrophila*.

COMPOSITION

Ingredients	Concentration
Sodium deoxycholate	50mg
Ampicillin	10mg

(per vial sufficient for 500 ml medium)

Applied in Medium










TM 1480 - RIPPEY CABELLI AGAR BASE

INSTRUCTION FOR USE

Rehydrate the contents of 1 vial aseptically with 5 ml of sterile distilled water. Aseptically add it to 500 ml of sterile, molten, cooled (45-50°C) TM 1480 - RIPPEY-CABELLI AGAR BASE. Mix well and pour into sterile petri plates.

STORAGE

Vials should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 Best Before	 QTY. Quantity	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 LOT/ B. NO. Lot / Batch Number	 Consults Instructions for Use	 QR Code	

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only
Revision: 3 March., 2022