

TS 216 – MEROPENEM SUPPLEMENT

INTENDED USE

For enrichment of Enterococci.

COMPOSITION

Ingredients	Concentration
Meropenem	1mg

(Per vial sufficient for 500 ml medium)

Applied in Medium

TM 1643 - VRE BROTH BASE (VANCOMYCIN RESISTANT ENTEROCOCCI BROTH)

TM 1642 - VRE AGAR BASE (VANCOMYCIN RESISTANT ENTEROCOCCI AGAR)

INSTRUCTION FOR USE

Rehydrate the contents of 1 vial with 5 ml of sterile distilled water and aseptically add to 1000 ml of sterile molten, cooled (45-50°C) TM 1642 - VRE AGAR BASE (VANCOMYCIN RESISTANT ENTEROCOCCI AGAR) along with 2 vials of Vancomycin Supplement or aseptically add 2 vials to 1000 ml of TM 1643 - VRE BROTH BASE (VANCOMYCIN RESISTANT ENTEROCOCCI BROTH). Mix well and dispense as desired.

STORAGE

Vials should be stored in sealed container at 2-8°C.

 Good Manufacturing Practices Certified	 For In Vitro Diagnostic Use	 Quantity	 Lot / Batch Number	 Catalogue Number	 Manufacturer
 Temperature Unit	 Authorized Representative MediStar GmbH Buckhorn 10, 48143 Münster, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

***For Lab Use Only**
Revision: 7 March., 2022