

TS 261 – β LACTAMASE I SUPPLEMENT (W/ 3300 IU OF PENICILLINASE)

INTENDED USE

For inactivation of penicillin group of antibiotics.

COMPOSITION

Ingredients	Concentration
Penicillinase	3300 IU of activity













(1 IU is defined as the amount of enzyme needed to hydrolyze 1 μ mole of Penicillin G (Penicillinase) per minute at 25°C and pH 7.0. 1 IU of Penicillinase corresponds to 600 Levy Units or 75 Pollock Units.)

INSTRUCTION FOR USE

Rehydrate the contents of 1 vial with 5 ml of distilled water. After reconstitution, filter sterilize the solution with 0.2 micron syringe filter. Add appropriate amount of solution depending on the application. Remaining sterilized solution can be stored at 2-8°C for 4 weeks.

STORAGE

Vials should be stored in sealed container at 2-8°C.

 Good Manufacturing Practices Certified	 For In Vitro Diagnostic Use	 Quantity	 Lot / Batch Number	 Catalogue Number	 Manufacturer
 Temperature Unit	 Authorized Representative MedNet GmbH Barkstrasse 10, 48163 Münster, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

***For Lab Use Only**
Revision: 16 June, 2023