

TS 265 – STERILE BETA LACTAMASE I SUPPLEMENT (READY TO USE)

INTENDED USE

For inactivation of penicillin.

COMPOSITION

Ingredients	Concentration
Penicillinase	Potency1 : 2000 Levy Units/ml/min Inactivation of Pen G2 :1,000,000 IU/ml/min










(per vial contains 20ml)

INSTRUCTION FOR USE

Sterile Beta lactamase I Supplement is an optimized ready to use solution that can be directly added to the test sample. The amount of product to be added to the test sample should be determined and set-up depending on the application, concentration of antibiotic to be inactivated and depending on the specific beta-lactam that has to be inactivated. Add appropriate amount of solution depending on the application. Remaining solution can be stored at 2-8°C. Aseptic techniques are to be followed throughout the procedure.

STORAGE

Vials should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 Best Before	 QTY. Quantity	 REF Cataloge Number	 Manufacturer
 Temperature Unit	 LOT/ B. NO. Lot / Batch Number	 Consults Instructions for Use	 QR Code	

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only
Revision: 4 March., 2022