

## TS 268 – MODIFIED CFC SELECTIVE SUPPLEMENT

### INTENDED USE

For selective isolation of *Pseudomonas* species.

### COMPOSITION

| Ingredients             | Concentration |
|-------------------------|---------------|
| Cetrimide               | 5mg           |
| Sodium fusidate         | 5mg           |
| Cephalothin sodium salt | 25mg          |

(Per vial sufficient for 500 ml medium)

### Applied in Medium

TM 2023 - CFC AGAR BASE (CEPHALOTHIN-SODIUM FUSIDATE-CETRIMIDE AGAR)

TM 2024 - CFC AGAR BASE (CEPHALOTHIN-SODIUM FUSIDATE-CETRIMIDE AGAR)

### INSTRUCTION FOR USE

Rehydrate the contents of one vial aseptically with 5 ml distilled water. Mix well and aseptically add it to 500 ml of sterile, molten, cooled (45-50°C), TM 2023 - CFC AGAR BASE (CEPHALOTHIN-SODIUM FUSIDATE-CETRIMIDE AGAR) / TM 2024 - CFC AGAR BASE (CEPHALOTHIN-SODIUM FUSIDATE-CETRIMIDE AGAR). Mix well and pour into sterile Petri plates.

### STORAGE

Vials should be stored in sealed container at 2-8°C.

|                                                                                                                                  |                                                                                                                                                                                            |                                                                                                            |                                                                                                           |                                                                                                                           |                                                                                                       |
|----------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------|
| <br>Good Manufacturing<br>Practices Certified | <br>For In Vitro Diagnostic Use                                                                         | <br>Quantity            | <br>Lot / Batch Number | <br>Catalogue Number                 | <br>Manufacturer |
| <br>Temperature Unit                          | <br>Authorized Representative<br><small>MedWet GmbH<br/>Südstraße 11<br/>48153 Münster, Germany</small> | <br>European Conformity | <br>QR Code            | <br>Consults Instructions<br>for Use | <br>Best Before  |

**NOTE:** Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

**\*For Lab Use Only**  
Revision: 7 March., 2022