

TMC 001B - PETRI PLATES CARRIER

INTENDED USE

Stainless steel petri plate dish carriers to facilitate the handling and the transport of petri plates around the laboratory and on site. It is useful for safe handling, transport and incubation of Petri dishes around the laboratory.

Product Features

The stainless steel construction, no corners to clean.

Swinging handle for easy manipulation is what the operators are asking for.

Used for sterilization and storage.

PACKING:

1X10 pack size

APPLICATION:

Useful for safe handling, transport and incubation of Petri dishes around the laboratory.

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only









