

TS 092 – LISTERIA SELECTIVE SUPPLEMENT-II

INTENDED USE

For the selective isolation and cultivation of *Listeria monocytogenes*.

COMPOSITION

Ingredients	Concentration
Acriflavin hydrochloride (Trypaflavin)	5mg
Nalidixic acid	20mg
Cycloheximide	25mg

(Per vial sufficient for 500 ml medium)

Applied in Medium

TM 1333 - BUFFERED LISTERIA ENRICHMENT BROTH BASE (BLE Broth Base) (as per FDA)

TM 1444 - LISTERIA SELECTIVE AGAR BASE

TM 761 - LISTERIA SELECTIVE BROTH BASE

TMV 761 - LISTERIA SELECTIVE BROTH BASE (VEG.)

INSTRUCTION FOR USE

Rehydrate the contents of one vial with 5 ml sterile distilled water. Mix well and aseptically add to 500 ml sterile, cooled (45-50°C) TM 1333 - BUFFERED LISTERIA ENRICHMENT BROTH BASE (BLE Broth Base) (as per FDA) / TM 761 - LISTERIA SELECTIVE BROTH BASE / TMV 761 - LISTERIA SELECTIVE BROTH BASE (VEG.) / TM 1444 - LISTERIA SELECTIVE AGAR BASE. Mix well and dispense as desired.

STORAGE

Vials should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 IVD For In Vitro Diagnostic Use	 QTY. Quantity	 LOT/ B. NO. Lot / Batch Number	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 EC REP Authorized Representative MedNet GmbH Buckhorn 10, 40143 Düsseldorf, Germany	 CE European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

***For Lab Use Only**
Revision: 7 March., 2022