

## TS 178 – ALBUMIN GLUCOSE SUPPLEMENT

### INTENDED USE

For rapid cultivation of *Mycobacterium tuberculosis*.

### COMPOSITION

| Ingredients        | Concentration |
|--------------------|---------------|
| Albumin fraction V | 1g            |
| Dextrose           | 1.500g        |
| Distilled water    | 20ml          |

(Per vial sufficient for 500 ml medium)

### APPLIED IN MEDIUM


TM 099 - DUBOS BROTH BASE

### INSTRUCTION FOR USE

Warm up the refrigerated supplement to room temperature. Aseptically add to sterile, cooled (45-50°C) 180 ml TM 099 - Dubos Broth Base. Mix well and pour into sterile screw capped tubes.

### STORAGE:

Vials should be stored in sealed container at 2-8°C.



**NOTE:** Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

**\*For Lab Use Only**  
Revision: 02 Feb., 2022