

TS 276 – V.P. SUPPLEMENT

INTENDED USE

For selective isolation & differentiation of *Legionella* species.

COMPOSITION

Ingredients	Concentration
Polymyxin B	50,000 IU
Vancomycin	1mg

(Per vial sufficient for 500 ml medium)

Applied in Medium


TM 983 - DIFFERENTIAL BUFFERED CHARCOAL YEAST EXTRACT AGAR BASE

INSTRUCTION FOR USE

Warm up the refrigerated solution to room temperature. Aseptically add the contents of 1 vial to 1000 ml of sterile, cooled (45-50°C) TM 983 - DIFFERENTIAL BUFFERED CHARCOAL YEAST EXTRACT AGAR BASE. Mix well and dispense into sterile tubes.

STORAGE

Vials should be stored in sealed container at 2-8°C.



NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

***For Lab Use Only**
Revision: 11 March., 2022